

CONSERVATIVE

ipoclor⁵

The excellent solutions for cleaning root canal

A choice based on scientific evidence IPOCLOR 5 Sodium hypochlorite at 5% active chlorine Ready to use solution based on sodium hypochlorite at 5% for the irrigation of root canals during reaming operations.

The high quality of glass packaging ensures greater stability of the product over time.

Characteristics

High capacity of dissolving the organic component because the high concentration (5.25) and the high pH;

Do not mix with acids because it would release toxic gases.

Do not inject - do not ingest Keep out of the reach of children.

Packaging

Multidose bottle 250 ml

Composition

Sodium hypochlorite at 5% active chlorine

Instructions for use

With a disposable syringe with blunt canal washing needle, aspirate the solution from the bottle and to wash canal.

CE

